



Bunch-by-bunch Beam Length Measurement Using Two-frequency System at SSRF

Hanjiao Chen, Liwu Duan, Yongbin Leng
SSRF BI Group





- Motivation
- Principles
- System Setup
- Beam Based Calibration
- Applications
- Next Works
- Summary

Motivation

For Users

- FB and Top-up make users happy usually
- Unstable beam conditions will damage data quality
 - Beam instabilities
 - The injection transient every 5mins (top-up)

Need better knowledge of instability phenomenon for optimization

For BI Engineers

- Only static information can be monitored during daily operation
- Injection transient is a very good window to learn more about machine

Get more information during routine injections

For Physicists

- A bunch is a basic unit for physics study
- Parameters of individual bunch is better than average values for physicists

Need parameters of every bunches

All requires




An eye on every bunches to *study* their behaviors



We are building a *6-dimensional bunch-by-bunch* diagnostic system

Motivation

Beam position (x,y)




Button pickups + BBB DAQ + $\Delta\Sigma$ algorithm

Resolution $\leq 10\mu\text{m}$

(Y. Yang, IPAC13; Y.B. Leng, IBIC13;
Y.B. Leng, IPAC15;)

Beam transverse size(horizontal , vertical)




SR light + PMT array + BBB DAQ
+ Gaussian fitting algorithm

(H. J. Chen, IPAC17; Under development)

6-dimensional bunch-by-bunch diagnostic system

Beam longitudinal phase (z)




Button pickups + BBB DAQ
+ Rising edge detection

Resolution $\leq 2\text{ps}$

(Y.M. Zhou, IPAC17; Poster TUPWC02 in IBIC17)

Beam longitudinal length



Button pickups + **Two-frequency system**

★ discussed in this talk

Streak Camera system at SSRF

Streak Camera Location:

- Hamamatsu C5680 Streak Camera
- Resolution: 2ps @high-speed sweep
50ps @low-Speed sweep
- Synchroscan rate: 125MHz (1/4 RF)
- Trigger rate: 2Hz (figure reputation rate)





Figure 1: Layout of the Optical Diagnostic Beamline (J. Chen, IPAC13)



individual bunch measurement long term measurement

Due to the limited pixels number

- Individual bunch measurement with very small time scale (sub us)
- Long term (ms) measurement with bunches averaging

Can't cover both

Streak Camera not Sufficient for 6D system

Our directions

- Bunch-by-bunch longitudinal oscillation studies
- One part of 6-dimensional bunch-by-bunch diagnostic system
- Intelligent trigger for beam abnormal status
- Bunch length control system for higher beam current in SSRF phase II, online monitor required

Our requirements

- **Real-time** bunch-by-bunch beam length measurement synchronized with other monitors (position, phase, etc)
- **High resolution** for individual bunch measurement (ps)
- Enough **measurement range** up to several thousands turns (ms)
- Easy to operate and connect to EPICS

Our limitations

- Off-line data analyze
- Low photons for bunch-by-bunch measurement
- Contradiction between enough data length and resolution of single bunch
- Complicated and stand-alone system,

Streak camera is good for dedicated machine study, not good for online monitoring. Frequency domain method (two-frequency) was adopted for 6D system. SC can be used to do calibration.

Principles of two-frequency method

- A typical bunch longitudinal charge distribution is **Gaussian**. Fourier component of ***mi-th harmonic*** is:

$$V(m_i \omega_0) = 2V_0 \exp\left(-\frac{m_i^2 \omega_0^2 \sigma_0^2}{2}\right) \quad (1)$$

RF frequency 499.654 MHz of SSRF

the bunch length
obtained with voltage ratio of two-frequency

$$\sigma = \sqrt{\frac{2}{m_2^2 \omega_0^2 - m_1^2 \omega_0^2} \ln\left(\frac{V_1(m_1 \omega_0)}{V_2(m_2 \omega_0)}\right)} \quad (2)$$

- V_1 and V_2 : the value at two *mi-th* harmonic frequencies
- $m_2 > m_1$

Basic Concept





Figure 3: Two-frequency method at KEK (Takao Ieiri, EPAC2000)
Two-frequency points: $1\text{GHz}(2f_{RF})$, $2.54\text{GHz}(5f_{RF})$
Successfully **average** bunch length measurement



'Old method'
Narrow band system,
slow DAQ

Figure 4: Two-frequency method at CLIC (H. H. Braun, EPAC1998)
Two-frequency points:
 $26\text{-}40\text{GHz}$; $60\text{-}90\text{GHz}$
Resolution around **0.7ps**
Single shot measurement

Bunch-by-bunch data acquisition techniques @SSRF




Figure 5: Hardware setup of the concept prototype of bunch-by-bunch processor (Yongbin Leng, IPAC2015)





Figure 6: Raw data and resampling points from oscilloscope (Yongbin Leng, IBIC2013)

Contribute to a
Bunch-by-bunch
two-frequency
bunch length monitor

Working Frequency Choice

The resolution:
$$\frac{\Delta\sigma}{\sigma} \simeq \frac{\sqrt{2}}{|m_1^2 - m_2^2| \omega_0^2 \sigma_0^2} \left[\left(\frac{\Delta V_1}{V_1} \right)^2 + \left(\frac{\Delta V_2}{V_2} \right)^2 \right] \quad (3)$$

➡ $m_1^2 - m_2^2 \uparrow$, resolution \downarrow



- Enough signal amplitude
- RF component limitation

Choose

$m_1=1$ and $m_2=6$

(about **500 MHz** and **3 GHz**)
as working frequencies

Figure 7: Beam Spectrum at storage ring of SSRF

Hardware Setup

HPF+LPF from WI company
Center Frequency: 500MHz
Bandwidth: ~ 500 MHz



Agilent DSO9064A Oscilloscope
Sampling rate: 5GS/s
Analog bandwidth: 600MHz







Figure 9: Data resampling with RF clock (499.654MHz)

LO signal as 5 times of RF signal
Generated from a programmable timing module EVAL-ADF4351




A mixer: m6-th signal to around 500 MHz with LO signal



Mini-Circuit
ZX05-42MH+

Figure 8: Block diagram of two-frequency bunch length measurement system
(L. W. Duan, *Nucl. Sci. Tech*, 2017)

Beam signal



Modification of Bunch Length Equation

$$\sigma = \sqrt{\frac{2}{m_2^2 \omega_0^2 - m_1^2 \omega_0^2} \ln \left(\frac{V_1(m_1 \omega_0)}{V_2(m_2 \omega_0)} \right)} \quad (2)$$

Formula modification due to two channels
transfer functions difference and *limited bandwidth*

$$\sigma = \sqrt{\frac{2}{m_2^2 \omega_0^2 - m_1^2 \omega_0^2} \ln \left(K_1 \frac{V_1}{V_2} + K_2 \right)} \quad (4)$$

No theoretical analytic solution for **K1** and **K2**
 The two coefficients can be calibrated by Steak Camera

m_2, m_1, ω_0 is theoretically knowable

$V1/V2$ is measured from two-frequency system

$$R_{12} = K_1 \frac{V_1}{V_2} + K_2 = \exp \left[\frac{\sigma^2 \cdot (m_2^2 \omega_0^2 - m_1^2 \omega_0^2)}{2} \right] \quad \text{Calculated from Steak Camera measurement}$$

Beam Based Calibration-Coefficient Calibration

Beam conditions

Single bunch filling
pattern

Bunch charge from
0.23nC to 6.05nC






Figure 10: Two-frequency system coefficient calibration by Steak Camera


Application: Bunch length measurement



Bunch length measurement with stored single bunch beam



Bunch length shortening during beam loss




**Bunch length has been measured
in very large dynamic range
(30pC to 6nC)**

**New method agreed with Streak
Camera very well**



Application: Longitudinal oscillation observation after Injections

Equivalent bunch length oscillation will be observed after injection
Macro bunch (refilled charge + stored charge)



- **Betatron oscillation** : From timing or amplitude mismatch of kickers
- **Synchrotron oscillation**: From phase or energy mismatch between injected charge and stored charge (double of the synchrotron frequency)


Application: Longitudinal oscillation observation after Injections




Equivalent bunch length oscillation after injection waveform can be calculated with the knowledge of

- ***Stored charge / Q_s***
- ***Refilled charge / Q_r***
- ***Synchrotron frequency / f_s***
- ***Longitudinal damping time / τ***
- ***Initial position in phase space of refilled charge / φ, z***


Or all above parameters can be retrieved by fitting the measured waveform




Application: Longitudinal oscillation observation after Injections




Filling pattern before and after injection



Measured bunch length of 3rd bunch trains




Bunch length variation of bunch #420 & # 430



FFT of the left waveforms

Next Work(1) – waveform fitting

- Find a way to do multi-parameters curve fitting
- Full information of refilled charge in phase space can be retrieved to guide injection optimization
- Injection quality can be evaluated precisely



Next Work(2) – DAQ upgrade

Random beam instability events can be captured with intelligent trigger




Figure 16: Two-frequency online system based on commercial DAQ

Next Work(3) – 6D BBB System integration





Figure 16: **6-dimensional bunch-by-bunch** diagnostic system architecture

Abnormal status ***capturing and analysis***
(Beam instabilities events, beam loss, etc.)

Next Work(4) – Multi-freq Schematic



Increase accuracy by adding more signal channels :
500MHz, 2GHz, 3GHz, 4.5GHz

Summary



- ***The two-frequency system*** with bunch-by-bunch capability is successfully implemented in SSRF.
- After the **beam based calibration using streak camera** as reference this system was used to measure bunch length in a very large charge dynamic range from 30pC to 6nC.
- This system also been used to monitor the **refilled bunch behavior during routine injection**. With further data analyzing the full information of refilled charge in phase space can be retrieved.
- **DAQ need to be upgraded** and **intelligent trigger mode** based on FPGA will be implemented. After upgrading this system will be more useful for physicist to capture bunch-by-bunch data when unstable beam condition shows up.

Acknowledge



- Appreciated for the support from National Natural Science Foundation of China (No. 11375255 and No. 11375254)
- Appreciated for the help from beam physics group of SSRF in design and analysis of beam experiment
- Appreciated for the help from beam operation group of SSRF in beam experiment



Thanks for your attention

Contact: chenhanjiao@sinap.ac.cn



中国科学院上海应用物理研究所
Shanghai Institute of Applied Physics, Chinese Academy of Sciences

Beam loss due to scrpation

