

BUNCH PATTERN VIA SINGLE PHOTON COUNTING AT SPEAR3

Abstract

A PicoHarp300 has been commissioned at SPEAR3 for Time-Correlated Single-Photon Counting. Initial results with a low dark-current PMT reveal a relatively large measurement error and charge spill on injection that need to be corrected for pump/probe operations. The PMT afterpulse can be removed by de-convolution. Measurements between adjacent topoff cycles clearly delineate charge injection into discrete buckets. A graphical interface in matlab extends the functionality of the PicoHarp300 for storage ring bunch pattern applications.

SLAC, SPEAR3 and Visible Beam Diagnostics



Aerial view of SLAC




SPEAR3 Light Source




Visible Diagnostic Beam Line




TCSPC Detector Configuration




Charge Spill and Afterpulse



Topoff Diagnostic and GUI Development



PicoHarp300 display interface



Summary

- PicoHarp300 TCSPC device at SPEAR3
- PMT with TTL conversion noisy
- Charge spill at injection identified
- afterpulse can be removed in software
- topoff measurements effective
- Matlab gui to extend capability